Postural Control in Bed

Recommended for the non-ambulant person with neurological impairment

Side Lying

stable side lying



Place folded sheet crosswise on bed. Bend hips and knees to 90°.




2 Roll person onto side. Support back with pillow tucked under crosswise sheet to hold in place.



3 Separate knees. Place roll between thighs, lower legs forward (in front) of roll.

> Loaded tissues can be relieved by alternating position side-to-side.



4 Place small cushion (pad) between feet to relieve pressure.



ISO 9001: 2008 COMPLIANT







Alternative Position

where abduction is a problem



1 Uncontrolled abducted lower limbs. Place folded sheet crosswise on bed.



2 Place legs over the abduction roll.



3 Adjust lateral pads according to individual. Place small cushion under lower legs to relieve pressure on heels. Support trunk laterally if necessary, by tucking a pillow under the crosswise sheet on each side.




Supine Lying

stable supine lying



1 Asymetrical uncontrolled (unsupported) lying. Place folded sheet crosswise on bed.



2 Align body posture. Place T Roll under knees with central roll flat.



3 Rotate central roll upwards and position in between thighs. Place small cushions or pillow under legs to raise heels of the bed.



 Support trunk laterally if necessary, by tucking a pillow under the crosswise sheet on each side. Knees must rest comfortably on roll.

C&S Seating Limited Phone 01424 853331 Fax 01424 854018 info@cands-seating.co.uk www.cands-seating.co.uk Company Reg 4718316

Offerfair Business Park **Ground Floor Front Office** 12-14 Maunsell Road St Leonards on Sea East Sussex TN38 9NN

Vat No 574 3844 14

REVISION AUGUST 2024 LB

